



# The ThermoSuit®



**ThermoSuit® in the "Lay Flat Mode"**



**Cross-Sectional View in the Inflated Mode**

The High-Technology Design of the ThermoSuit® Lower Fluid Delivery System Incorporates:

- Head positioning indicator
- Porous lining for uniform fluid distribution
- Velcro sealing for easy routing of patient lines in and out
- Expandable gussets to accommodate a wide variety of patients



**Complete ThermoSuit® System**

**No Hospital Should be Without The ThermoSuit® System!**



170 Kinnelon Road • Kinnelon, NJ 07405

Tel: (973) 283-2800 • Fax: (973) 283-2910 • Email: [sales@life-recovery.com](mailto:sales@life-recovery.com) • Web: [life-recovery.com](http://life-recovery.com)

PN 51008, Ver A